

Collagen Membrane

Collagen Membrane

- Easy manipulation
- Dual-sided usage
- Barrier function lasting for 4-6 months

*Not all products listed in this catalog are available in all countries.

www.dentium.com

CAT-M0502(Rev.1, 1101)

Distributor **Dentium** For Dentists By Dentists 3105 Trade Tower 159, Samsung-dong, Gangnam-gu, Seoul, South Korea 135-729 T +82-2-501-8560 / F +82-2-567-9578


Manufacturer **GENOSS** 1F Gyeonggi R&DB Center 906-5 Iui-dong, Yeongtong-gu, Suwon-si, Gyeonggi-do, Korea T +82-31-888-5100 / F +82-31-888-5105

Dentium
For Dentists By Dentists

Collagen Membrane

RESORBABLE MEMBRANE



Application of Collagen Membrane



- Thinner membrane (300 μ m) with multiple layers for easy manipulation and sufficient mechanical strength in surgery.

- Resorption period of 4-6 months to provide enough time for stabilizing graft materials and supporting bone growth.

- Multiple-layered structure enables more effective bone regeneration by sparing enough space for hard tissue formation and facilitates proliferation of osteoblast.




SEM

Product


Product	Size	Thickness
GCM1020	10X20mm	0.3mm
GCM1520	15X20mm	
GCM2030	20X30mm	
GCM3040	30X40mm	

Preclinical Data

- Rabbit Calvaria Model, 6-12 weeks



- Degradation character in collagenase solution



Clinical Case

